

Leading Edge Genomic Services & Solutions for FHS

Date : 25 Feb 2021

Elaina Lo



澳門大學
UNIVERSIDADE DE MACAU
UNIVERSITY OF MACAU



健康科學學院
Faculdade de Ciências da Saúde
Faculty of Health Sciences

A low-angle photograph of a modern, multi-story building with a light-colored, textured facade and numerous windows. The building is partially obscured by large, diagonal blue and white graphic elements that sweep across the frame from the bottom left towards the top right. The Novogene logo is visible on the lower part of the building.

Novogene
诺禾致源

Content

1. **Service Introduction**
 - Comprehensive NGS Services
2. **Services Workflow**
 - Extensive Quality Control and Standard Workflow
3. **Overview of Novogene**
 - Providing Advanced Genomic Solutions

A low-angle photograph of a modern, multi-story building with a light-colored, metallic-looking facade and large windows. The building is set against a clear blue sky. Several thick, diagonal blue and light blue lines are overlaid on the image, extending from the top left towards the bottom right.


Novogene
诺禾致源

Content

1. **Service Introduction**
 - Comprehensive NGS Services
2. Services Workflow
 - Extensive Quality Control and Standard Workflow
3. Overview of Novogene
 - Providing Advanced Genomic Solutions

Comprehensive Research Services and Solutions

Using our array of state-of-the-art **Illumina**, **Pacbio** and **Nanopore** systems, Novogene provides comprehensive genomics, transcriptomics and proteomics research solutions.




First intelligent multi-product NGS delivery platform

Falcon system


Launched in 2020

Fully automated production line → reduce manual errors and basis

Safe and accurate one-stop ground-breaking solution



DNA Sequencing



Sequencing Region:
Whole Genome

Sequencing Depth:
>30X-50X

Identify all variants: SNPs,
INDELs, SV, CNV

Sequencing Region:
Whole Exome

Sequencing Depth:
>100X-200X

Identify most variants: SNPs,
INDELs in **coding region**, CNV
for paired cancer samples

Sequencing Region:
Specific Region (customized)

Sequencing Depth:
>500X

Identify most variants: SNPs,
INDELs in **specific region**

Comparison of WGS, WES & TRS


Sequencing strategy	Applications	Advantages	Disadvantages
WGS	<ul style="list-style-type: none"> Genetic disease study Cancer research Human population evolution DNA biomarkers Pharmacogenomics Evolutionary and demographic history Biodiversity 	<ul style="list-style-type: none"> No need of prior knowledge of disease Focus on non-coding area mutations unbiased Detection of all variant types including CNV & SV Joint Analysis of Multi-Groups 	<ul style="list-style-type: none"> High cost Data interpretation is difficult Low frequency mutations cannot be detected
WES	<ul style="list-style-type: none"> Genetic disease study Cancer research Human population evolution Focus only on coding areas for mutation 	<ul style="list-style-type: none"> Low cost Good detection sensitivity 	<ul style="list-style-type: none"> Can't detect non-coding region and viral SV detection is incomplete and inaccurate
TRS	<ul style="list-style-type: none"> Cancer research Human population studies Linkage analysis for inherited diseases Discovery of biomarkers and therapeutic targets Candidate genes screened for pre-sequencing in large sample sizes High-depth validation of mutation information 	<ul style="list-style-type: none"> Low cost Can detect low frequency mutations Offers more targeted detection of mutations 	<ul style="list-style-type: none"> Needs knowledge of gene of interest Probes need to be customized and arrive with long delivery times

Whole Genome Sequencing (WGS)

Service Parameter

Platform	Illumina
Read length	2×150 bp
Turnaround time	About 15 working days
Standard analysis	additional 10 working days
Advanced analysis	upon request

Project Workflow




Whole Genome Sequencing (WGS)

Sample Requirements

Sample Type	Amount (Qubit®)		Volume	Concentration	Purity (NanoDrop™/Agarose Gel)
	Strongly Recommended	Required			
Genomic DNA	≥0.4μg	≥ 0.2μg	≥ 20μL	≥ 10ng/μL	OD260/280 = 1.8-2.0, no degradation, no contamination
Genomic DNA (PCR-free)	≥ 3μg	≥ 1.5μg	≥ 15μL	≥ 10ng/μL	
FFPE	≥ 1.6μg	≥ 0.8μg	≥ 15μL	≥ 10ng/μL	Fragment ≥ 1500 bp


Whole Genome Sequencing (WGS)

Library construction



Whole Genome Sequencing (WGS)

Standard Bioinformatics



Whole Genome Sequencing (WGS)

Advanced analysis contents


Single gene disorder	Cancer	Complex disease
<ol style="list-style-type: none">1. Filter screened mutations in dbSNP and 1000 Genome database;2. Screen the variation of exomic or splicing site;3. Amino acid conservation prediction by SIFT and Polyphen;4. Gene function and pathway analysis;5. Encode/miRbase annotation of non-coding region.	<ol style="list-style-type: none">1. Oncogene/ Antioncogene/Susceptibility gene screening;2. Cancer tissue purity and ploidy analysis;3. Mutation characteristics analysis;4. Tumor driver gene analysis;5. Hypermutation analysis;6. Pathway enrichment analysis;7. Encode/miRbase annotation of non-coding region.	<ol style="list-style-type: none">1. Filter screened mutations in dbSNP and 1000 Genome database;2. Screen the variation of exomic or splicing site;3. Amino acid conservation prediction by SIFT and Polyphen;4. Gene function and pathway analysis;5. Encode/miRbase annotation of non-coding region;6. Novo mutation screening7. Analysis of gene significance;8. Protein-protein interaction analysis

Whole Exome Sequencing (WES)

Service Parameter

Platform	Illumina
Exome Capture	Agilent SureSelect Kit V6
Read length	2×150 bp
Turnaround time	about 18 working Days
Standard analysis	additional 4 working days
Advanced analysis	upon request

Project Workflow




Whole Exome Sequencing (WES)

Sample Requirements (Agilent V6)

Sample Type	Amount (Qubit®)		Volume	Concentration	Purity (NanoDrop™/Agarose Gel)
	Strongly Recommended	Required			
Genomic DNA	≥0.8μg	≥ 0.4μg	≥ 15μL	≥ 20ng/μL	OD260/280 = 1.8-2.0, no degradation, no contamination
FFPE	≥ 1.6μg	≥ 0.8μg	≥ 15μL	≥ 20ng/μL	Fragment ≥ 1000 bp
cfDNA/ctDNA	≥ 80ng	≥ 40ng	≥ 15μL	≥ 0.5ng/μL	Fragment in multiples of 170 bp, no genomic contamination


Whole Exome Sequencing (WES)

Library construction



Whole Exome Sequencing (WES)

Standard Bioinformatics




Whole Exome Sequencing (WES)

Advanced analysis contents

Single gene disorder	Cancer	Complex disease
<ol style="list-style-type: none">1. Filter screened mutations in dbSNP and 1000 Genome database;2. Screen the variation of exomic or splicing site;3. Amino acid conservation prediction by SIFT and Polyphen;4. Gene function and pathway analysis.	<ol style="list-style-type: none">1. Oncogene/ Antioncogene/Susceptibility gene screening;2. Cancer tissue purity and ploidy analysis;3. Mutation characteristics analysis;4. Tumor driver gene analysis;5. Hypermutation analysis;6. Pathway enrichment analysis.	<ol style="list-style-type: none">1. Filter screened mutations in dbSNP and 1000 Genome database;2. Screen the variation of exomic or splicing site;3. Amino acid conservation prediction by SIFT and Polyphen;4. Gene function and pathway analysis;5. Novo mutation screening6. Analysis of gene significance;7. Protein-protein interaction analysis

RNA Sequencing



RNA Family	Classification	Name	Percentage
	Coding RNA	Messenger RNA, mRNA	2%-5%
	Noncoding RNA	Long noncoding RNA, lncRNA	~5%
		Small RNA	
		microRNA, miRNA	
		Small interference RNA, siRNA	
		Piwi-protein interacting RNA, piRNA	80%-85%
		Ribosomal RNA, rRNA	
	Transfer RNA, tRNA		~15%


The transcriptome is the set of all RNA molecules in one cell or a population of cells. Based on next generation sequencing technique, transcriptome sequencing is capable of fully and quickly acquiring the overall transcripts information in the specific tissues or organs of eukaryotic species under a certain status.

Comparison of mRNA Seq, small RNA Seq & lncRNA

Sequencing strategy	Applications	Remarks
mRNA Seq	<ul style="list-style-type: none">• Pathological mechanism• Tumor-subtypes classification• Molecular markers• Human evolution• Drug target• Clinical diagnostics• Biology Development	<ul style="list-style-type: none">• Captures both known and novel features included gene expression structure and quantification• Identify biomarkers across the broadest range of transcripts• Can detect alternative Splicing
Small RNA Seq (size < 200nt) e.g., miRNA, piRNA	<ul style="list-style-type: none">• Expression quantification of small RNA transcripts• Function verification, such as gene knockout, over-expression of miRNA genes• Advanced Analysis: miRNA target gene verification, piRNA identification and expression quantification	<ul style="list-style-type: none">• Understand how post-transcriptional regulation contributes to phenotype• Identify novel biomarkers• Only available for eukaryotic species with reference genome
lncRNA Seq (size ≥ 200nt) e.g., lincRNA, NAT	<ul style="list-style-type: none">• Expression Quantification of lncRNA transcripts• Gene or RNA subcellular Localization and Expression• Function Verification, such as gene knockout, over-expression of lncRNA genes• Protein Interaction	<ul style="list-style-type: none">• Included mRNA information

Experimental Planning of RNA Seq

1. Number and type of replicates
 - 3 biological replicates being the **minimum** for any inferential analysis
 - More replicate is better, especially for clinical samples
 - Try to avoid pooling of individuals/experiments
2. Avoiding confounding
 - Ensure samples in each condition are all the same sex, age, litter, and batch, if possible.
3. Addressing batch effects
 - Avoid batch or split replicates across batches
4. Using same sample for sequencing and validation




mRNA Sequencing

Service Parameter

Platform	Illumina
Read length	2×150 bp
Turnaround time	15-18 working days
Standard analysis	additional 5-15 working days
Advanced analysis	upon request

Project Workflow




mRNA Sequencing

Sample Requirements

Library Type	Sample Type	Amount (Qubit®)		RIN	Volume	Concentration	Purity (NanoDrop™/Agarose Gel)
		Strongly Recommended	Required				
Eukaryotic (PolyA enrichment)	Total RNA (Eukaryotic)	≥0.8μg	≥ 0.4μg	≥ 6.8 (animals)/ 6.3 (plants & fungi)	≥ 20μL	≥ 20ng/μL	OD260/280 ≥ 2.0, no degradation, no contamination
Eukaryotic (Directional)	Total RNA (Eukaryotic)	≥1.6μg	≥ 0.8μg				
Prokaryotic (Directional, rRNA removal)	Total RNA (Prokaryotic)	≥6μg	≥ 3μg	≥ 6	≥ 20μL	≥ 50ng/μL	


mRNA Sequencing

Library construction



mRNA Sequencing

Standard Bioinformatics




Small RNA Sequencing

Service Parameter

Platform	Illumina
Read length	SE50
Library type	18-40 bp insert sRNA library
Turnaround time	20-22 working days
Standard analysis	additional 5-10 working days
Advanced analysis	upon request

Project Workflow




Small RNA Sequencing

Sample Requirements

Sample Type	Amount (Qubit®)		RIN	Volume	Concentration	Purity (NanoDrop™/Agarose Gel)
	Strongly Recommended	Required				
Total RNA	≥4μg	≥ 2μg	≥ 7.5 (animals)/ 7 (plants & fungi)	≥ 20μL	≥ 50ng/μL	OD260/280 ≥ 2.0, OD260/230 ≥ 2.0, no degradation, no contamination
Total RNA (Exosome)	≥40ng	≥ 20ng	Peak Range: 25-200 nt, FU≥ 10, with no peak >2000nt			

Small RNA Sequencing

Standard Bioinformatics



Whole Genome Bisulfite Sequencing

Service Parameter

Platform	Illumina
Read length	PE 150 bp
Turnaround time	25 working days (< 10 samples)
Standard analysis	additional 10 working days

Project Workflow




Whole Genome Bisulfite Sequencing

Sample Requirements

Sample Type	Amount (Qubit®)		Volume	Concentration	Purity (NanoDrop™/Agarose Gel)
	Strongly Recommended	Required			
Genomic DNA	≥5μg	≥ 2.5μg	≥ 15μL	≥ 20ng/μL	OD260/230 = 0-3, no degradation, no color, no RNA or protein contamination
Low-input Genomic DNA	≥ 0.4μg	≥ 0.2μg	-	≥ 20ng/μL	


Whole Genome Bisulfite Sequencing

Library construction



Whole Genome Bisulfite Sequencing

Standard Bioinformatics




A low-angle photograph of a modern, multi-story building with a light-colored, metallic-looking facade and large windows. The building is set against a clear blue sky. The Novogene logo is visible on the lower part of the building.

Novogene
诺禾致源

Content

1. Service Introduction
 - Comprehensive NGS Services
2. **Services Workflow**
 - Extensive Quality Control and Standard Workflow
3. Overview of Novogene
 - Providing Advanced Genomic Solutions

Sample Shipping Workflow (Macau)



Sample Shipping Workflow (Macau)

Dry ice provider (Taobao):

- 1) 源太干冰
- 2) 火星干冰

Please order 1 day in advance and ship to ZhuHai

Shipping address of Tianjin:

Novogene BuildingB07,Venture Headquarter Base,Wuqing Development Area,Tianjin,China.


天津市武清区创业总部基地B07栋

Contact: Sample Receiving Department


联系人：天津科服收样组

Tel: +86 18522699037


Extensive Quality Control



Sample QC




Library QC



We have strict quality control steps, only qualified library can be used for sequencing and analysis. The quality control include volume, concentration and library peak shape.

Data Quality





Content

1. Service Introduction
 - Comprehensive NGS Services
2. Services Workflow
 - Extensive Quality Control and Standard Workflow
3. Overview of Novogene
 - Providing Advanced Genomic Solutions

A low-angle photograph of a modern, multi-story building with a light-colored, metallic-looking facade and numerous windows. The building is partially obscured by large, diagonal blue and white graphic elements that sweep across the frame from the bottom left towards the top right. The sky is a clear, bright blue.

Novogene Mission



Founded in March 2011

With cutting-edge genomic technologies and high-performance computing, Novogene is dedicated to improving life science research and human healthcare.

Novogene
诺禾致源

Executive Team




Dr. Ruiqiang Li, Ph.D.

Founder

Dr. Ruiqiang Li, is widely recognized as one of the world's leading experts in genomics and bioinformatics, owning 17 patents and has authored over 100 publications that have been cited more than 42,000 times. Best known for developing the software SOAP (Short Oligonucleotide Analysis Package) for ultra-fast sequence mapping, variation detection, and de novo genome assembly, Dr. Li has focused the company on being a full-service provider, performing NGS services and providing bioinformatic analysis to fully empower and enable our customers to move forward rapidly and effectively in their scientific endeavors.

Novogene Strategy



Novogene Global Presence



Largest Range of Platforms to Serve Diverse Needs

>280,000 Human WGS/Year



NovaSeq6000



Pacbio Sequel
& Sequel II



Nanopore



HiSeq X



HiSeq
2000/2500/4000



Orbitrap
Exploris™ 480



DA8600



Q Exactive™ HF-X



S5XL



NextSeq
500



Miseq



QTRAP® 6500+

First intelligent multi-product NGS delivery platform

Falcon system

Launched in 2020

Fully automated production line → **reduce manual errors and basis**

Safe and accurate one-stop ground-breaking solution





Super Computing and Data Center


280,000 Human WGS Analysis / Year
1,000 Human WGS Analysis / Day

High Performance Computing

Capability	Capacity
Memory Size	410 TB
Computing power	1727 T flops
Storage	60.2 PB



Novogene
诺禾致源



Summary of Novogene Competitive Advantages

- ✓ **Corporation** – Global company, multi-geographic laboratories with multinational talents.
- ✓ **Platform** – The largest sequencing center equipped with the latest and comprehensive platforms.
- ✓ **Team** – 1,800 professionals, dedicated and specialised groups.
- ✓ **Quality** – Clinical standards (CAP, CLIA & ISO13485), automatic workflow, LIMS.
- ✓ **Experience** – 28,000+ projects and one million samples in 2017, 26,000+ global customers.
- ✓ **Innovation** – Sustainable development powered by R&D (over £15 Million), comprehensive product portfolio
- ✓ **Research** – 580+ Publications, global collaboration

Thank You !

www.novogene.com

Contact: Elaina Lo

elaina.lo@novogene.com

asia_hmt@novogene.com

